

High Throughput Screening tools for water quality monitoring

B. van Vugt^a, E. Simon^a, M. Naderman^a, J. Legradi^b, M. Lamoree^b, J. Legler^b, B. van der Burg^a
^aBioDetection Systems bv, Amsterdam; ^bInstitute for Environmental Studies (IVM), Vrije Universiteit Amsterdam

Introduction


For water quality monitoring, chemical analysis is widely used and well-grounded in regulatory frameworks. However, the practical use of chemical analysis is severely hampered by the need for *a priori* selection of compounds.

As an alternative, several effect-based bioanalytical tools for water quality evaluation have been developed in FS8.1.1 and 8.1.4. Effect-based assays have the advantage that they provide information on the biological effect of compounds, irrespective of their chemical nature.

Aim

Recently, two case-studies [1] have been performed where a representative set of water samples was analysed using a broad panel of CALUX reporter gene assays (figure 1).

The aim of the current study was to demonstrate the added value of such an effect-based analysis. To this end, 81 prioritized compounds were screened on the CALUX panel, and the resulting activity profile was compared to the CALUX profile of the water samples (figure 2).




Approach

The water samples and a set of 81 routinely monitored substances (figure 3) were measured on the human reporter gene CALUX panel (figure 4), thus generating a toxicity profile.

PRIORITY HAZARDOUS	PRIORITY	UNDER REVIEW	OTHERS
1 Anthracene	1 Abotolol	1 AMPA	Pharmaceuticals, list of compounds provided by Swiss and Italian water companies
2 Benzo(a)pyrene	2 Atrazine	2 Bentazon	
3 Benzo(b)fluoranthene	3 Benzene	3 Bisphenol-A	
4 Benzo(g,h,i)perylene	4 Chlorfeniphos	4 Dicofol	
5 Benzo(k)fluoranthene	5 Chlorpyrifos-ethyl	5 EDTA	
6 Indeno(1,2,3-cd)pyrene	6 1,2-Dichloroethane	6 Free cyanide	
7 C10-13-chloroalkanes	7 Dichloromethane	7 Glyphosate	
8 Cadmium and its compounds	8 Di(2-ethylhexyl)phthalate (DEHP)	8 Mecoprop	
9 cadmium chloride	9 Fluoranthene	9 Musk xylene	
10 Endosulfan	10 Heptachlor	10 PCB118	
11 Hexachlorobenzene	11 Heptachlor epoxide	10 PCB126	
12 Hexachlorobutadiene	12 Lead and its compounds	10 PCB128	
13 Hexachlorocyclohexane	13 lead chloride	10 PCB156	
14 Mercury and its compounds	14 Nickel and its compounds	11 PFOS	
15 methylmercury(II) chloride	15 nickel (II) chloride	12 Quinoxifen	
16 mercuric chloride	16 Nickel (II) chloride	13 TCDD	
17 Nonylphenol technical mixture	17 Pentachlorophenol		
18 Pentachlorobenzene	17 Simazine		
19 PBDE 100	18 Trichlorobenzenes		
20 PBDE 47	19 Trichloromethane + chloroform		
21 Tributyltin-cation / hydride	20 Trifluralin		

Figure 3. Selected compounds: monitored or under review



Screening results


The compounds and the samples activated similar molecular pathways: induction of xenobiotic metabolism, endocrine pathways, and reactive mode-of-action. However, the particular assays activated within these categories were different (table 1). For example, the water samples consistently showed activity on pathways related to oxidative stress (Nrf2) and lipid metabolism (PPAR), while these pathways were hardly activated by the routinely monitored compounds.

Molecular pathways	Activated by monitored compounds:	Activated by Water Samples case studies:
Xenobiotic Metabolism:	DR CALUX PXR CALUX	DR CALUX PXR CALUX PPAR CALUX
Endocrine (nuclear hormone receptors):	ERa CALUX antiAR CALUX antiPR CALUX	ERa CALUX antiAR CALUX GR CALUX
Reactive MoA:	p53 CALUX	Nrf2 CALUX
Cell Viability:	Cytotox CALUX	Cytotox CALUX

Table 1. Activated molecular pathways and corresponding CALUX assays for both datasets

Activity profiles

Table 2 shows CALUX panel profiles of the compounds (left) and samples (right). The black boxes indicate matches between the two datasets (ERa, antiAR, PXR, DR, Cytotox), while the yellow boxes indicate mismatches (PPARs, Nrf2).



Conclusions

The CALUX profile of a set of water samples was compared to the profile of 81 priority pollutants. Although there was considerable overlap between the activated pathways in both datasets, some pathways activated by the water samples were clearly absent in the monitored compounds dataset. This suggests that these compounds cannot fully account for all activities observed in water samples.

The majority of the 81 selected compounds showed activity in at least one CALUX assay. A select panel of CALUX assays could be used to replace laborious chemical analysis; this would enable the detection of the priority pollutants, as well as other compounds that might be present in water samples. The present study clearly shows the added value of effect-based monitoring compared to standard chemical analysis.

References

[1] Escher BI et al, Environ Sci Technol, 48,1940-56. 2014.



BioDetection Systems bv
 Email: info@bds.nl
 Internet: www.bds.nl

