



BioDetection Systems

10th Biodetectors Meeting 2017 Sorrento, Italy

State of the Art Biodetectors Special emphasis on health monitoring


**Abraham Brouwer,
CEO of BDS & MLS, Amsterdam
Professor of Environmental Toxicology & Ecogenomics,
VU University Amsterdam
Managing Director of BE-Basic Consortium on biobased economy**



Core activity: Develop and apply novel bioassays for health & safety assessment




Safety assessment Authenticity testing Screening for bioactive compounds



Projects, partnerships

Demonstration studies



Novel bioanalytical solutions

Safety assessment:

- ❖ Food,
- ❖ Feed,
- ❖ Water,
- ❖ Health
- ❖ Environment
- ❖ Chemicals
- ❖ Consumer products

CALUX bioassays make use of cellular signal transduction pathways i.e., AOP-based bioassays


Dioxins, PAH

Toxic metals

PPAR activators

Endocrine Disruptors

Genotoxic compounds



Dioxin receptor

ROS, Nrf pathway



PPAR receptors

Hormone receptors

P53 pathway

Key Benefits:

- ❖ High predictivity of health related-effects
- ❖ Good estimate of total effect from mixtures
- ❖ Can predict unknown effects of chemicals
- ❖ Can discover unknown chemicals in matrices
- ❖ Level of precision similar to instrumental methods
- ❖ Low cost, high capacity, easy to operate




Best Usage:

- ❖ Most valuable tool for (human) biomonitoring
- ❖ Powerful screening tool for safety assessment e.g food, water
- ❖ Good *in vitro* alternative for chemical safety assessment



Nuclear receptors			Signaling pathways			Controls		
name	status	cell	name	status	cell	name	status	cell
DR CALUX	✓	U2OS					✓	U2OS
PAH CALUX	✓						✓	all
ER CALUX	✓						✓	all
ERalpha CALUX	✓						✓	all
ERbeta CALUX								
ERalpha CALUX								
ERbeta CALUX								
AR CALUX								
PR CALUX								
GR CALUX								
TR CALUX								
RAR CALUX								
PPARγ1 CALUX								
PPARγ2 CALUX								
PPARα CALUX	✓	U2OS	STAT CALUX	✓	U2OS			
PPARδ CALUX								
LXR CALUX								
PXR CALUX								
VDR CALUX								
MR CALUX	✓	U2OS						

- Acute toxicity
- Oxidative stress
- AhR pathway
- Endocrine effects/EDCs
- obesogens
- Reproductive effects
- Genotoxicity/carcinogenicity
- Metabolism
- etc



CALUX: n=28
Agonist/antagonist: 25x2=56 assays


Using the Ah-Receptor as AOP in DR-CALUX for monitoring dioxin toxicity by mixtures of dioxin-like compounds





Comparison of Dioxin analysis in e.g food DR CALUX[®] vs HR-GC/MS

Total dioxin-levels (PCDDs, PCDFs and dioxin-like (dl)-PCBs) in fishoil





BDS develops bio-based detection methods and applies those in a wide range of sectors

Environment
technology




Human
epidemiology



Medicin



Technique and
processes




Quality
Management



Feed



Food





Pharma



AOP Example 2: Hormone Receptors and endocrine disruptors

Using adverse endocrine/reproductive pathways to develop a set of ER,AR,PR,GR,TR, etc-CALUX systems in a human cell line for monitoring endocrine disruptors and reprotoxic chemicals





Predictability of CALUX for Reproductive hazard (risk) identification

	COMPOUND	TOXICITY	EST diff	zebrafish		CALUX panel	CALUX with PBPK	cyp19	PREDICTION
1	Cyclosporin A (CSA)	developmental (immuno) toxicant	differentiation effect	no effect	reproductive	anti AR, weak antiPR and GR and ESRE		no effect	positive
2	Monoethylhexylphthalate (MEHP)	male reproductive organ malformations	differentiation effect	developmental toxicant	neurodevelopmental	PPARg and PPARalpha agonist		no effect	positive
3	Sodium valproate (VPA)	neurodevelopmental toxicant	differentiation effect	developmental toxicant	neurodevelopmental	weakly positive in many assays, consistent with HDAC inhibition		no effect	positive
4	D-mannitol (DML)	negative control	no effect	no effect	neurodevelopmental	negative		no effect	negative
5	Flusilazole (FLU)	craniofacial and axial skeletal malformations	differentiation effect	developmental toxicant	reproductive	cytotoxic antiPR/antiGR weak DR/PAH		inhibitor at high conc in H295R	positive
6	Glufosinate ammonium (GPA)	neurodevelopmental toxicant	no effect	no effect	neurodevelopmental	negative		no effect	negative
7	Methoxyacetic acid (MAA)	growth and developmental retardation	differentiation effect	developmental toxicant	neurodevelopmental	negative		no effect	positive
8	Retinoic acid (RA)	neural crest cell migration affected	differentiation effect	developmental toxicant	neurodevelopmental	strong RAR/RXR activity		inhibitor in H295R	positive
9	Diocetyl tin dichloride/ dichlorodioctylstannane(DO TC)	developmental (immuno)toxicant	cytotoxic	no effect	neurodevelopmental	cytotoxic, antiprogesterin, stress-related pathways		inhibitor in H295R	positive
10	Endosulfan (ESF)	neurotoxicant	cytotoxic	developmental toxicant	neurodevelopmental	cytotoxic, ER, antiAR, antiGR		inhibitor in H295R	positive
11	Diethylstilbestrol (DES)	transplacental carcinogen	cytotoxic	developmental toxicant	neurodevelopmental	strong estrogen: antiAR, antiPR, stress- and genotoxicity		no effect	positive
12	Methylmercury chloride (MMC)	neurodevelopmental toxicant	cytotoxic	developmental toxicant	neurodevelopmental	stress-related pathways affected, estrogen, GR agonist		Inducer in H295R, inhibitor in HPMS	positive

Piersma et al. 2013 Reprod Toxicol. 38:53-64.



Surface water quality



Waste water treatment



Human monitoring



Safety & quality of Food packaging materials




Anabolic steroid abuse


AOP Example 3: p53-Calux and carcinogens

Using a cancer-based adverse outcome pathway to develop a p53-CALUX in a human cell line for chemical carcinogens and mutagens testing

P53 "the guardian of the genome":



Typical p53 CALUX responses by chemical carcinogens




	p53 CALUX (+/-S9)
Sensitivity (%)	82
Specificity (%)	90

***Validated using the ECVAM recommended list of 61 compounds
(Kirkland et al 2008. Mutat. Res. 653, 99–108)***

AOP Example 4: PPAR-Calux and obesity & diabetes-type 2

Using an metabolic syndrome/obesogen-based adverse outcome pathway to develop a set of PPAR-CALUX systems in a human cell lines for testing of chemicals that can induce obesity and diabetes-type 2





CALUX panel can measure AOP-based signatures

Option: identify 21st century priority chemicals

no activity
 EC10 = 1E-3M
 EC10 = 1E-7M

compound	Cytox10%	Cytox50%	ERA	ERA-anti	ERB	ERB-anti	AR	AR-anti	PR	PR-anti	GR	GR-anti	TRb	RAR	LXR	PPARa	PPARg	DR	PAH	Hif1a	TCF	AP1	ESRE	NFKB	Nf2	p21	p53	
Chlordane	-5.5	-5	-6.9					-6.5		-6.5		-6																
DDT	-4.5	-4.2	-6.5		-5.8			-7		-6		-5.5				-5						-4.7	-3.5					
Dieldrin	-3		-5.8					-7		-7		-5				-5.5												
Endrin			-5.5					-7		-7																		
Heptachlor	-5	-4.5	-7.2					-7		-6																		
Hexachlorobenzene			-6.5					-6		-6																		
Mirex																												
Toxaphene	-5	-4.8	-5.5		-5.5			-6.5		-6.5		-5.5																
PCB118	-4.5							-7		-6.5																		
PCB126	-4.8	-4.4						-6.5		-6																		
PCB128								-7		-6.5																		
PCB156	-4.5	-4	-6					-6		-6																		
TCDD																												
Furan																												
dibenzo[a,h]anthracene	-4				-7.5																							
dibenzo[a,h]pyrene			-7																									
benzo[a]pyrene			-6		-3.9			-6.5		-6																		
tributyltinacetate								-7																				
methylmercury(II)chloride	-7	-6.5																										
	-5.8	-5.6			-6.4							-6.0																
Lead chloride	-3.5	-3																										
Mercuric chloride	-4.8	-4.8																										
Cadmium chloride	-4.9	-4.7																										
Cobaltous chloride	-3.9	-3.4																										
Copper chloride																												
copper sulfate	-3.4	-3.2																										
Zinc sulphate	-4.3	-4.1																										
Sodium arsenite	-5.4	-5.2																										
Nickel(II)chloride	-3.5	-3																										
chromium(vi)oxide	-5	-4.7																										

Dirty Dozen POPs: endocrine, dioxin receptor



Additional POPs: dioxin receptor, stress pathways

Heavy metals: acute toxicity, stress pathways

BDS: Bioactivity: discovery of highly potent anti-tumor, Antibiotic, anti-oxidant and anti-obesity activities in plants & biomass



anti-AR CALUX



BioDetection Systems

BDS is supported by grants from BE Basic






Applications of CALUX bioassays in health monitoring (focus on metabolic syndrome)

Some examples of BDS involvement in health monitoring

- **Association between EDCs and equine metabolic syndrome**
- **(collaboration with Dr Mickelson and Dr. McCue of University of Minnesota College of Veterinary Medicine)**
- **Newborns and Genotoxic exposure risks in EU-Newgeneris**
- **(collaboration with many parties within the EU-Newgeneris consortium)**

Health monitoring: Endocrine disruptors and Metabolic syndrome development

- Dioxins cause wasting syndrome, *AhR* maybe involved in metabolic syndrome development
- Endocrine *ER*, *TR* and *PPAR* pathways are involved in regulation of fat and energy metabolism and insuline sensitivity
- Dioxins interfere in *ER* and *TR* pathways;
- Dioxins and endocrine disruptors suspected of supporting Diabetes type II development



The association between endocrine disrupting chemicals and equine metabolic syndrome

S.A. Durward-Akhurst¹, E.M. Norton¹, N. Schultz¹,
R. Geor², J. Mickelson¹, M.E. McCue¹

¹University of Minnesota College of Veterinary Medicine, St Paul, MN

²University of Massey, Turitea, NZ


10th Biodetectors Sorrento 2017



UNIVERSITY OF MINNESOTA
Driven to DiscoverSM

Equine Metabolic Syndrome

- Increased adipose deposition
 - Regionally
 - Generalized obesity
- Insulin dysregulation
- Predisposition to laminitis



10th Biodetectors Sorrento 2017



UNIVERSITY OF MINNESOTA
Driven to DiscoverSM

Hypothesis



1. Endocrine disrupting chemicals are associated with the EMS phenotype
2. The *AHR* and/or *ER* genotype of an individual modulates the metabolic response secondary to exposure to EDCs

10th Biodetectors Sorrento 2017

We ALL are Affected! Get the Facts!

Hide Your Kids... Hide Your Wife... Cause Koppers Superfund is KILLING EVERYONE Up In HERE!

Study design

- 161 Morgan horses
 - 18 farms
- 140 Welsh ponies
 - 14 farms

- Phenotypic measurements:

- Fasting glucose and insulin
- Post OST glucose and insulin
- Triglycerides
- ACTH
- Leptin
- Adiponectin
- NEFAs



Statistics

- Non-normally distributed data was transformed
- Univariate linear model
- Response variables:
 - Fasting glucose and insulin
 - Post OST glucose and insulin
 - Triglycerides
 - ACTH
 - Leptin
 - Adiponectin
 - NEFAs


10th Biodefactors Sorrento 2017



UNIVERSITY OF MINNESOTA

Driven to DiscoverSM

AHR-EDC results




Median: 19.29 pg/g fat
Mean: 38.51 pg/g fat
Range: 0.59-536.36 pg/g fat
Below LOD: 131 samples

n = 301



ER-EDC results



Median: 5.50 pg/ml

Mean: 101.91 pg/ml

Range: 4.35-15000.00 pg/ml

Below LOD: 119 samples

n = 276



Results

- signifies $p > 0.05$

EMS phenotype	AHR-EDC	ER-EDC
Resting insulin	-	$p = 0.003$
INS-OST	-	$p = 0.002$
Resting glucose	$p = 0.042$	$p = 0.002$
GLU-OST	-	$p = 0.012$
NEFA	$p = 0.047$	-
Triglycerides	$p = 0.011$	-
Adiponectin	-	-
Leptin	-	-
ACTH	-	-




Conclusions

- EDCs acting through the ER and AHR are associated with EMS in Welsh Ponies and Morgans
- EDCs likely explain some of the unexplained environmental variance in EMS phenotype



UNIVERSITY OF MINNESOTA

Driven to DiscoverSM



A recent study has also found evidence for a possible association between Dioxins and Human Metabolic Syndrome



Perinatal exposure to dioxins and dioxin-like compounds and infant growth and body mass index at seven years: A pooled analysis of three European birth cohorts

Nina Iszatt [a](#), Hein Stigum [a](#), Eva Govarts [b](#), Lubica Palkovicova Murinova [c](#), Greet Schoeters [b,d,e](#), Tomas Trnovec [c](#), Juliette Legler [f,1](#), Cathrine Thomsen [g](#), Gudrun Koppen [b](#), Merete Eggesbø [a,□](#)

RESULTS

At 7 years, dioxins exposure was associated with a statistically significant increase in BMI in girls (adjusted estimate for BMI units $\beta=0.49$, 95% CI: 0.07, 0.91) but not in boys ($\beta=-0.03$, 95% CI: -0.55, 0.49) (p-interaction=0.044). Furthermore, girls had a 54% (-6%, 151%) increased risk of overweight at 7 years (p-interaction = 0.023).

Conclusions

Perinatal exposure to dioxin and dioxin-like compounds was associated with increased early infant growth, and increased BMI in school age girls. Studies in larger sample sizes are required to confirm these sexspecific effects.

Please cite this article as: Iszatt, N., et al., Perinatal exposure to dioxins and dioxin-like compounds and infant growth and body mass index at seven years: A pooled analysis of..., Environ Int (2016), <http://dx.doi.org/10.1016/j.envint.2016.04.040>



Applications of CALUX bioassays in health monitoring (early life stage exposure & effects)

Some examples of BDS involvement in health monitoring

- **Association between EDCs and equine metabolic syndrome**
- **(collaboration with Dr Mickelson and Dr. McCue of University of Minnesota College of Veterinary Medicine)**
- **Newborns and Genotoxic exposure risks in EU-Newgeneris**
- **(collaboration with many parties within the EU-Newgeneris consortium)**

Hypothesis to be tested:

Maternal exposure to dietary compounds with carcinogenic and immunotoxic properties results in *in utero* exposure and molecular events in the unborn child leading to increased risk of cancer and immune disorders later in childhood.

Existing mother-child cohorts will be used while new bio banks will be set-up

Overall goal:


Development and application of two categories of biomarkers in relation to dietary exposure and childhood disease.



- 1 - biomarkers of exposure to chemicals with carcinogenic and immunotoxic properties
- 2 - biomarkers of pre-carcinogenic and immunotoxic effects



Box Plot of DR CALUX results on blood plasma from Newgeneris Cohorts





BioDetection Systems

Outcome of NewGeneris Study



- ❖ Pathway specific bioassays are valuable for human monitoring
- ❖ Small volume sample analysis of human plasma is feasible with CALUX bioassays
- ❖ 11 papers published:
The NewGeneris human early lifestage epidemiology studies show associations between exposure to dioxins and/or EDCs (especially with cord serum) and adverse Health outcome in children, in particular:
 - **Associations between DR-CALUX responses and childhood leukemia**
 - Associations between DR-CALUX responses and low birth weight; and shorter gestational age
 - Associations between DR-CALUX responses and changes in AGD in young boys
 - Prenatal exposure to DR-CALUX responses via food is associated with effects on the immune system functions at 1 and 3 year old children

NewGeneris results:

Possible relation between dioxin exposure and incidence of childhood leukemia

Global Gene Expression Analysis in Cord Blood Reveals Gender-Specific Differences in Response to Carcinogenic Exposure In Utero

Kevin Hochstenbach¹, Danitsja M. van Leeuwen¹, Hans Gmuender⁴, Ralf W. Gottschalk¹, Martinus Løvik⁵, Berit Granum⁵, Unni Nygaard⁵, Ellen Namork⁵, Micheline Kirsch-Volders⁶, Ilse Decordier⁶, Kim Vande Loock⁶, Harrie Besselink², Margareta Törnqvist⁷, Hans von Stedingk⁷, Per Rydberg⁷, Jos C.S. Kleinjans¹, Henk van Loveren^{1,3}, and Joost H.M. van Delft¹ *Cancer Epidemiol Biomarkers Prev* 2012;21:1756-1767. Published OnlineFirst August 9, 2012.

Methods: Global gene expression was applied in umbilical cord blood samples, the **CALUX-assay** was used for measuring dioxin(-like), androgen(-like), and estrogen(-like) internal exposure, and acrylamide–hemoglobin adduct levels were determined by mass spectrometry adduct-FIRE-procedureTM. To link gene expression to an established phenotypic biomarker of cancer risk, micronuclei frequencies were investigated

Conclusions/Impact: This study reveals different transcriptomic responses to environmental carcinogens between the sexes. In particular, **male-specific TNF-alpha-NF-kB signaling upon dioxin exposure and activation of the Wnt-pathway in boys upon acrylamide exposure might represent possible mechanistic explanations for gender specificity in the incidence of childhood leukemia**



NewGeneris results:

Env. carcinogens, endocrine disruptors may mechanistically contribute to carcinogen-induced childhood leukemia

Micronuclei in Cord Blood Lymphocytes and Associations with Biomarkers of Exposure to Carcinogens and Hormonally Active Factors, Gene Polymorphisms, and Gene Expression: The NewGeneris Cohort

Domenico Franco Merlo,¹ Silvia Agramunt,² Lívia Anna,³ Harrie Besselink,⁴ Maria Botsivali,⁵ et al., Environ Health Perspect 122:193–200, february 2014; ; <http://dx.doi.org/10.1289/ehp.1206324>

CALUX-relevant part of Results

Gene expression levels were significantly lower for 11 genes in association with the highest versus lowest category of plasma AR CALUX® (chemically activated luciferase expression for androgens) (8 genes), ER α CALUX® (for estrogens) (2 genes), and DR CALUX® (for dioxins).

Conclusion: We measured *in utero* exposure to selected environmental carcinogens and circulating hormonally acting factors and detected associations with Micronuclei frequency in newborns circulating T lymphocytes. The results highlight mechanisms that may contribute to carcinogen-induced leukemia and require further research.





We are happy to discuss any options for future collaboration



Thank you for your attention!