

# Applying the CALUX bioassay panel: GREEN toxicology

**BioDetectors Conference 2018**


BARBARA VAN VUGT-LUSSENBURG,  
BIODETECTION SYSTEMS




This project has received funding from the Bio Based Industries Joint Undertaking under the European Union's Horizon2020 research and innovation programme under agreement No 745450.

- Effect-based risk assessment, focusing on molecular initiation events
  - chemical – receptor interaction (agonism/antagonism)
  - chemical – cell signaling pathway interaction (activation/inhibition)
- Suitable for single compounds or complex mixtures
- >30 CALUX assays available

Cell line	Endpoint
ERa CALUX (ago/anta)	Estrogen receptor (ant)agonists
AR CALUX (ago/anta)	Androgen receptor (ant)agonists
PR CALUX (ago/anta)	Progesterone receptor (ant)agonists
GR CALUX (ago/anta)	Glucocorticoid receptor (ant)agonists
TRb CALUX (ago/anta)	Thyroid receptor (ant)agonists
RAR CALUX	Retinoic acid receptor agonists
LXR CALUX	Liver X receptor agonists
PXR CALUX	Pregnane X receptor agonists
PPARa CALUX	Peroxisome proliferator activated receptor agonists
PPARg2 CALUX	Peroxisome proliferator activated receptor agonists
PPARd CALUX	Peroxisome proliferator activated receptor agonists
AhR CALUX	Aryl Hydrocarbon receptor agonists
Hif1a CALUX	Chemical hypoxia response
TCF CALUX	wnt/TCF pathway activation
AP-1 CALUX	AP1 pathway activation / cell cycle control
ESRE CALUX	Endoplasmic reticulum stress
NFkB CALUX	Activation of NF-kB pathway (immune response)
Nrf2 CALUX	Oxidative stress
p21 CALUX	Transcription of p21 inhibitor of cell cycle progression
p53 CALUX	p53-dependent pathway activation / genotoxicity
Cytotox CALUX	Cytotoxicity



- CALUX panel analysis results in a hazard profile:







# Increasing the throughput: automation

To increase throughput AND improve quality: automate the most time consuming steps

- cell seeding (cell dispenser)
- compound dilution (liquid handling robot)
- compound exposure (liquid handling robot)
- luminescence detection (luminometer + stacker)



CELL SEEDING



EXPOSURE




DETECTION


## Manual vs robot

Identical EC50-, SD- and RLU values:


Dose-response curve DR CALUX®



Dose-response curve AR CALUX®



Same 'hits' are found by robot as by two experienced technicians:





- ↓ Labour time
- ↓ Interindividual variation
- ↓ Human errors
- ↑ Reproducibility
- ↑ Data quality

# Increasing the throughput: miniaturization

96-wells -> 384-wells plates resulted in:


- 2 -> 6 compounds per plate
- 8 -> 16 concentrations per compound
- 1 control -> full dose-response curve per plate
- 100 -> 200 dose-response curves per day
  - broader concentration range
  - more controls
  - more reliable and accurate data
  - automated analysis possible
  - lower sample volumes





- Toluene and N-methyl-2-pyrrolidone (NMP) are two of the most toxic solvents currently in use.
- The primary aim of the ReSolve project is to replace these solvents with safe, biobased alternatives derived from non-food carbohydrates.

## The concept





# Establishing an integrated testing strategy to evaluate safety issues

- *in vitro*: human CALUX cell-based reporter gene assay panel
  - major types of toxicity
  - rapid
- *in silico*: DIAMONDS, CoMSAS, TTC
  - combining existing data: read-across
  - structure/activity evaluation
- Analyse solvents currently in use. Their profile can be seen as a point of departure: toxicity profiles of novel candidates should be more favourable.

# Modify CALUX protocols to enable screening of volatile solvents

- Expose solvents at 1% v/v: exposure concentration  $\sim 0.1\text{M}$  instead of  $1\text{E}^{-4}\text{M}$  (standard procedure)
- Cover the CALUX assay plates with non-breathing seals after exposure to prevent solvent evaporation
- Use  $\text{CO}_2$ -independent cell culture medium to allow cell survival in a closed environment (due to the seal)



# CALUX panel results for reference solvents




BioDetection Systems

toxicity

xenobiotics  
sensor

stress pathways


compound	toxicity										xenobiotics sensor				stress pathways										
	Cytotox20%	ERa-anti	AR	AR-anti	PR	PR-anti	GR	GR-anti	TRb	TRb-anti	RAR	PXR	PPARa	PPARd	PPARg	Ahr	Hif1a	TCF	Ap1	ESRE	NFKB	Nrf2	p21	p53 GENTOX +/-S9	
NMP	-1.4					-1.7		-1.5						-2.2		-2.0		-2.5	-1.5					-2.0	
Toluene	-2.2											-1.0											-2.0		
Xylene	-3.0											-3.1													
ether 1,4-dioxane						-0.9						-2.1													
acetonitrile	-0.7											-1.2													
2-(ethoxymethyl)tetrahydrofuran	-1.6					-2.0						-2.1													
1,2-dichlorobenzene	-3.8																								
propylene carbonate	-1.5																								
Gamma-valerolactone (GVL)	-1.4											-2.0													
DMSO																									
sulfolane																									
1,3-Dimethyl-2-imidazolidinone	-1.0									-2.5			-2.5			-2.5			-1.5					-2.0	
1,3-Dimethyl-3,4,5,6-tetrahydro-2(1H)-pyrimidinone	-1.1									-2.5			-2.6			-2.6		-2.6						-2.0	
DMF																									-1.3
DMAC													-2.2			-2.2		-2.5	-1.7						-2.0
benzene												-1.7													



# CALUX panel results for solvent candidates



- Some candidates look promising
- Some are (cyto)toxic and/or active on too many assays
- Some possess undesirable properties: endocrine activity





- The CALUX panel was successfully optimized for this specific class of compounds (volatile solvents)
- The CALUX results were used to identify candidates with a favourable safety profile, for additional (tox) studies and structural optimization
- Some of the results obtained so far had not been predicted based on structural information and physicochemical properties of the candidates
- The ReSolve project is therefore a good example of the added value of **effect-based** bioassays
- Currently, the CALUX *in vitro* data is being combined with *in vivo* / *in silico* data to identify missing endpoints in the current panel

# Acknowledgements

- WFBR

David Franciolus

Linda Gootjes

Shanmugam Thiyagarajan

Frits van der Klis

Daan van Es

- UOY

Fergal Byrne

James Sherwood

Thomas Farmer

- NOVA

Asta Partanen

Lara Dammer

Ángel Puente

- TNO

Dinant Kroese

Lisette Krul (left TNO)

Harrie Buist

Marie Meima

- BDS

Barbara van Vugt

Bart van der Burg

Hai-Yen Man

Faried Niamut

Leon Schweisshelm