


Comprehensive safety screening of substances in European waste and surface waterbodies using a large panel of 25 CALUX bioassays

Harrie Besselink¹, Barbara van Vugt¹, M. Naderman¹, E. Simon², Bram Brouwer^{1, 3}, Bart van der Burg¹
¹BioDetection Systems bv, Amsterdam, the Netherlands; ³Oekotoxzentrum/Centre Ecotox, Zürich, Switzerland;
³Vrije Universiteit Amsterdam, Amsterdam, the Netherlands

Introduction

- Surface and (bio-based) waste waters might contain various unknown toxicants, metabolites, degradation products, which are not included in current regulations.
- Direct measurement of possible effects of chemical contaminants in extracts of water samples using bioassays is expected to deliver a better approach for water quality assessment compared to measurement of a limited number of individual chemicals.
- To demonstrate bioassays as complementary tools to chemical analytics we aimed at a comprehensive toxicity screening of the currently monitored substances in European waterbodies and compared the identified toxic endpoints to case studies on effect-based water quality screening.

Approach




- Toxicity profile of the routinely/occasionally monitored substances in European surface waters (n=66)
- Identification of the toxic pathways picked up by these substances
- Comparison of the toxicity pathways induced by the routinely monitored substances with relevant toxicity pathways for water quality assessment (suggested by various case studies [1-3])

Materials and Methods


In vitro reporter gene assay panel (25 CALUX) was first compiled covering a wide range of toxic endpoints, that are assumed to be indicative for water samples.

Figure 1 – The CALUX assay principle



The assays have been automated using a compact liquid handling system. Dilution series of 16 concentrations in triplicate are tested for each chemical to derive PC₁₀ (agonism; 10% induction relative to reference) and PC₈₀ (antagonism; 80% induction relative to reference) values.

Figure 2 – Workflow of automated compound screening



References

- van der Linden S. Applicability of functional genomics tools for water quality assessment. 2014. PhD Thesis, VU University, Amsterdam
- Leusch FD et al., 2014. Water Res. 2014. 49:300-15.
- Escher BI et al., 2014. Environ Sci Technol 48(3):1940-56.

Results

Substance	CALUX Assays																								
	DR CALUX	PXR CALUX	PPARα CALUX	PPARγ CALUX	ERα CALUX	ERβ CALUX	AR CALUX	PR CALUX	GR CALUX	Nr1h2 CALUX	PC10 CALUX	PC80 CALUX	Cytotoxic CALUX	PC10 CALUX	PC80 CALUX	PC10 CALUX	PC80 CALUX	PC10 CALUX	PC80 CALUX	PC10 CALUX	PC80 CALUX	PC10 CALUX	PC80 CALUX		
1. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
2. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
3. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
4. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
5. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
6. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
7. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
8. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
9. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
10. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
11. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
12. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
13. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
14. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
15. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
16. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
17. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
18. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
19. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
20. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
21. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
22. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
23. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
24. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	
25. Atrazine	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	

Table 1: Toxicity screening of the set of substances routinely/occasionally monitored in European water bodies. Values represent logarithmic PC₁₀ (agonistic assays) and PC₈₀ (antagonistic assays) concentrations.

Conclusions

- Compound screening revealed the importance of ENDOCRINE - (particularly the activation of the ERα-, anti-AR, anti-PR receptors), GENOTOXICITY (p53-CALUX) and XENOBIOTIC METABOLISM-related (PXR-CALUX) pathways of the routinely monitored substances (Table 1).
- This profile of pathway effects is consistent with only part of the profile observed in water samples (Figure 3) indicating that not all substances are picked up by the compound panel of the monitoring programs.

Activated by monitoring compounds	Activated by water samples [1,2,3]
DR CALUX	✓
PXR CALUX	✓
ERα CALUX	✓
Anti-AR CALUX	✓
Anti-PR CALUX	✓
GR CALUX	✓
P53 (-S9) CALUX	✓
Nr1h2 CALUX	✓
Cytotoxic CALUX	✓

Figure 3 – Identified toxic pathways of the chemically monitored substances and compared with the pathways indicated by Dutch and Australian case studies [1-3] on effect-based screening of water samples

This study demonstrates that effect-based methods could complement conventional chemical analysis in water quality monitoring as prescreening techniques by

- identifying toxic "hotspots" for further investigation,
- assessing the effect of the entire mixture of compounds present in waters and therefore,
- reduce uncertainty in safety evaluation.